

Remember –
Optimize delivery of
enteral nutrition –
Elevate head of bed to 45°*

45° Optimal elevation

30° Minimum elevation

*Unless contraindicated

CANADIAN GUIDELINES FOR
Achieving Best Practice
CRITICAL CARE NUTRITION

A collaboration of  Critical Care
Nutrition and  Nestlé
Nutrition

311-05/09E

Remember –
Optimize delivery of
enteral nutrition –
Elevate head of bed to 45°*

45° Optimal elevation

30° Minimum elevation

*Unless contraindicated

CANADIAN GUIDELINES FOR
Achieving Best Practice
CRITICAL CARE NUTRITION

A collaboration of  Critical Care
Nutrition and  Nestlé
Nutrition

311-05/09E

Remember –
Optimize delivery of
enteral nutrition –
Elevate head of bed to 45°*

45° Optimal elevation

30° Minimum elevation

*Unless contraindicated

CANADIAN GUIDELINES FOR
Achieving Best Practice
CRITICAL CARE NUTRITION

A collaboration of  Critical Care
Nutrition and  Nestlé
Nutrition

311-05/09E

Remember –
Optimize delivery of
enteral nutrition –
Elevate head of bed to 45°*

45° Optimal elevation

30° Minimum elevation

*Unless contraindicated

CANADIAN GUIDELINES FOR
Achieving Best Practice
CRITICAL CARE NUTRITION

A collaboration of  Critical Care
Nutrition and  Nestlé
Nutrition

311-05/09E

Remember –
Optimize delivery of
enteral nutrition –
Elevate head of bed to 45°*

45° Optimal elevation

30° Minimum elevation

*Unless contraindicated

CANADIAN GUIDELINES FOR
Achieving Best Practice
CRITICAL CARE NUTRITION

A collaboration of  Critical Care
Nutrition and  Nestlé
Nutrition

311-05/09E

Remember –
Optimize delivery of
enteral nutrition –
Elevate head of bed to 45°*

45° Optimal elevation

30° Minimum elevation

*Unless contraindicated

CANADIAN GUIDELINES FOR
Achieving Best Practice
CRITICAL CARE NUTRITION

A collaboration of  Critical Care
Nutrition and  Nestlé
Nutrition

311-05/09E